



Virant Diagnostics, Inc.

11002 Veirs Mill Rd, Suite 404
Wheaton, MD 20902
Phone: (877) 888-2973, Fax: (888) 713-3456
CAP #: 954036801, CLIA #: 21D2184276
www.virantdx.com

ImmunoCAP IgE and IgG4 Allergy Test Requisition Form

Place Barcode Label Here

Blood Specimen Information

Collection Date: ___/___/___
Collection Time: ___:___ AM PM

PATIENT INFORMATION

Last Name:	First Name:	Date of Birth:
Email:	Phone:	Gender: <input type="checkbox"/> M <input type="checkbox"/> F
Address:	City: State:	Zip:

DIAGNOSIS CODES

- | | | |
|---|---|--|
| <input type="checkbox"/> T78.1XXA: Adverse food reactions | <input type="checkbox"/> T78.00XA: Anaphylactic reaction due to unspecified food, initial encounter | <input type="checkbox"/> L50.0: Allergic urticaria |
| <input type="checkbox"/> J30.1: Allergic rhinitis due to season | <input type="checkbox"/> J30.81: Allergic rhinitis due to animal | <input type="checkbox"/> H10.45: Chronic allergic conjunctivitis |
| <input type="checkbox"/> Z88.0: Allergy to penicillin | <input type="checkbox"/> Z91.030: Bee venom allergy | <input type="checkbox"/> Z91.038: Other insect venom allergy |
| <input type="checkbox"/> Z91.010: Allergy to peanuts | <input type="checkbox"/> Z91.011: Allergy to milk products | <input type="checkbox"/> Z91.012: Allergy to eggs |
| <input type="checkbox"/> Z91.018: Allergy to other foods | <input type="checkbox"/> E55.9: Vitamin D deficiency | <input type="checkbox"/> Other: |

INSURANCE AND PAYMENT INFORMATION

Bill Insurance (Attach copy of insurance card, front and back) Bill Client Bill Patient (Cash/Check/Credit Card) Other:

Primary Plan Name:	Policy Holder Name:
Policy #:	Group #:
Secondary Plan Name:	Policy Holder Name:
Policy #:	Group #:

HEALTHCARE PROVIDER INFORMATION

Provider Name:	NPI:
Organization:	Phone:
Address:	City: State: Zip:

Provider Authorization Signature: _____ Date: ___/___/___

SAMPLE REQUIREMENTS

Specimen Collection and Storage: For each panel ordered, completely fill one red-top tube or gel barrier tube. Separate serum by centrifugation and transfer the serum into a plastic transport tube. Store at 2°C to 8°C up to 7 days, or else at -20°C. Avoid repeated freezing and thawing.

Specimen Transport: Ship specimens overnight at room temperature (15°C to 25°C). If above room temperature, then the package must include an ice pack.

TOTAL IgE (82785)

Total IgE 82785 x 1

sIgE ALLERGEN PANELS WITH REFLEX (86003)

If whole allergen tests positive (>0.10 kU/L) – reflexes to the component panels. Please see the next page for detailed allergens included in each panel. All panels can be modified according to individual preference.

<input type="checkbox"/> Region 2 Respiratory Panel (86003 x40)	<input type="checkbox"/> Common Food Panel (86003 x 12)
<input type="checkbox"/> Nut Panel (86003 x9)	<input type="checkbox"/> Early Childhood Profile (<2 year old) (86003 x16)
<input type="checkbox"/> Fish Panel (86003 x7)	<input type="checkbox"/> Shellfish Panel (86003 x7)
<input type="checkbox"/> Venom Panel (86003 x5)	<input type="checkbox"/> Legume Panel (86003 x5)
<input type="checkbox"/> Occupation Panel (86003 x12)	<input type="checkbox"/> Drug Panel (86003 x11)
<input type="checkbox"/> Custom Panel	<input type="checkbox"/> Custom Panel
<input type="checkbox"/> Custom Panel	<input type="checkbox"/> Custom Panel

INDIVIDUAL ALLERGEN sIgE TESTS

Over 650 allergens are available for specific IgE and IgG4 testing. Please check our website for a comprehensive list of allergen tests available. To view the full catalogue, please visit www.virantdx.com/testing-solutions/immunocap-allergy-testing/.

<input type="checkbox"/> Egg white	<input type="checkbox"/> Milk	<input type="checkbox"/> Peanut	<input type="checkbox"/> Latex	<input type="checkbox"/> Wheat	<input type="checkbox"/> Gluten
<input type="checkbox"/> Soybean	<input type="checkbox"/> Milk, boiled	<input type="checkbox"/> Sesame	<input type="checkbox"/> Shrimp	<input type="checkbox"/> Lobster	<input type="checkbox"/> Crab
<input type="checkbox"/> Other:	<input type="checkbox"/> Other:	<input type="checkbox"/> Other:	<input type="checkbox"/> Other:	<input type="checkbox"/> Other:	<input type="checkbox"/> Other:

sIgE COMPONENT TESTS (86008)

Performs the allergen component tests without the whole allergen. Over 100 component tests are available for specific IgE and IgG4 testing.

<input type="checkbox"/> sIgE Dog Component Panel 86008 x 6	<input type="checkbox"/> sIgE Cat Component Panel 86008 x 4
<input type="checkbox"/> sIgE D. pteronyssinus Component Panel 86008 x 4	<input type="checkbox"/> sIgE Hazelnut Component Panel 86008 x 4
<input type="checkbox"/> sIgE Cashew Component Panel 86008 x 1	<input type="checkbox"/> sIgE Brazil nut Component Panel 86008 x 1
<input type="checkbox"/> sIgE Walnut Component Panel 86008 x 2	<input type="checkbox"/> sIgE Peanut Component Panel 86008 x 8
<input type="checkbox"/> sIgE Egg Component Panel 86008 x 2	<input type="checkbox"/> sIgE Milk Component Panel 86008 x 3
<input type="checkbox"/> sIgE Soy Component Panel 86008 x 3	<input type="checkbox"/> sIgE Wheat Component Panel 86008 x 4
<input type="checkbox"/> sIgE Sesame Component Panel 86008 x 1	<input type="checkbox"/> sIgE Gliadin Component Panel 86008 x 1
<input type="checkbox"/> sIgE Venom Component Panel 86008 x 8	<input type="checkbox"/> sIgE Shrimp Component Panel 86008 x 1
<input type="checkbox"/> sIgE Alternaria alternata Component Panel 86008 x 1	<input type="checkbox"/> sIgE Common silver birch Component Panel 86008 x 4
<input type="checkbox"/> sIgE Timothy grass Component Panel 86008 x 6	<input type="checkbox"/> sIgE Giant ragweed Component Panel 86008 x 1
<input type="checkbox"/> sIgE Alpa-gal Panel 86008 x 1	
<input type="checkbox"/> Other:	<input type="checkbox"/> Other:
<input type="checkbox"/> Other:	<input type="checkbox"/> Other:

IgG4 TEST MENU (86001)

Specific IgG4 tests are not FDA cleared. They are laboratory derived tests and have been validated by Virant Diagnostics for clinical purposes. Most allergens or components can be tested for both sIgE and sIgG4.

<input type="checkbox"/> IgG4 Timothy grass	<input type="checkbox"/> IgG4 Hazelnut	<input type="checkbox"/> IgG4 Sesame seed	<input type="checkbox"/> IgG4 Peanut
<input type="checkbox"/> IgG4 Giant ragweed	<input type="checkbox"/> IgG4 Cashew	<input type="checkbox"/> IgG4 Milk	<input type="checkbox"/> IgG4 rAra h 1
<input type="checkbox"/> IgG4 Common silver birch	<input type="checkbox"/> IgG4 Pistachio	<input type="checkbox"/> IgG4 Egg white	<input type="checkbox"/> IgG4 rAra h 2
<input type="checkbox"/> IgG4 Oak	<input type="checkbox"/> IgG4 Walnut	<input type="checkbox"/> IgG4 Cat dander	<input type="checkbox"/> IgG4 rAra h 3
<input type="checkbox"/> IgG4 D. pteronyssinus	<input type="checkbox"/> IgG4 Pecan	<input type="checkbox"/> IgG4 Dog dander	<input type="checkbox"/> IgG4 rAra h 6
<input type="checkbox"/> IgG4 D. farinae	<input type="checkbox"/> IgG4 Brazil nut	<input type="checkbox"/> Other:	<input type="checkbox"/> Other:

ALLERGY PANELS

Region 2 Respiratory Panel	D. farinae, D. pteronyssinus, Cat dander, Dog dander, Mouse urine proteins, German Cockroach; (Tree Pollens x 10) Box-elder, Common silver birch, Cottonwood, Elm, Maple leaf sycamore, Mountain juniper, Mulberry, Oak, Walnut, White ash; (Grass pollens x 8) Bahia grass, Bermuda grass, Cocksfoot (Orchard), Johnson grass, Meadow fescue, Meadow grass (Kentucky blue), Rye-grass, Timothy grass; (Weed pollens x 10) Cocklebur, Common pigweed, Common ragweed, Giant ragweed, Goosefoot (Lamb's quarters), Mugwort, Nettle, Plantain (English) Ribwort, Rough marshelder, Sheep sorrel, (Mold x 6) Alternaria alternata, Aspergillus fumigatus, Cladosporium herbarum, Mucor racemosus, Penicillium chrysogenum, Stemphylium herbarum
Early Childhood Allergy Profile	D. farinae, D. pteronyssinus, Cat dander, Dog dander, Mouse urine proteins, German cockroach, Cladosporium herbarum, Alternaria alternata, Egg white, Peanut, Soybean, Milk, Shrimp, Walnut, Cod, Wheat
Common Food	Egg white, Cow's milk, Peanut, Soybean, Sesame, Wheat, Walnut, Cashew nut, Shrimp, Cod, Salmon, Tuna
Nut Panel	Almond, Brazil nut, Cashew nut, Hazel nut, Peanut, Pecan nut, Pine nut, Pistachio, Walnut
Fish Panel	Cod, Octopus, Pacific squid, Plaice, Salmon, Trout, Tuna
Shellfish Panel	Blue mussel, Clam, Crab, Lobster, Oyster, Scallop, Shrimp
Legume Panel	Chickpea, Lentil, Pea, Soybean, White bean
Venom Panel	Common wasp, Honey bee, Paper wasp, White-faced hornet, Yellow hornet
Drug Panel	Amoxicilloyl, Ampicilloyl, Cefaclor, Chlorhexidine, Gelatin bovine, Insulin human, Penicilloyl G, Penicilloyl V, Pholcodine, Morphine, Suxamethonium
Occupation Panel	Bougainvillea, Cotton seed, Ethylene oxide, Ficus, Formaldehyde/Formalin, Isocyanate HDI, Isocyanate MDI, Isocyanate TDI, Ispaghula, Phthalic anhydride, Silk, Sunflower seed

Custom designed panels and other regional respiratory panels are available upon request.

FOR LABORATORY USE ONLY

Accession #:	Patient ID:
Client/Physician ID:	Date Received: ___/___/___
	Time Received: _____: _____ <input type="checkbox"/> AM <input type="checkbox"/> PM

Please visit our webpage at www.virantdx.com/testing-solutions/immunocap-allergy-testing/ for more information.
Contact us at (877) 888-2973 or allergylab@virantdx.com for any inquiries.